(D)

TECHNOBASE® 6000/8000

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(GB)

TECHNOBASE® 6000/8000

Universal resins for veterinary medicine

User instructions

Material

TECHNOBASE 6000 and TECHNOBASE 8000 are rapid-curing resins on a methyl methacrylate basis supplied as powder and liquid.

Indication

Claw treatment:

Immobilisation of claws during inflammatory diseases. Fixation of claw tips.

Treatment of bone fractures in large and small animals:

Extra-cutaneous osteosynthesis splinting after Becker, transversal fixation by means of Kirschner drill wires, radius and ulna fractures, tibia-fibula fractures, pedal bone fractures, strengthening of dressings.

Properties

TECHNOBASE 6000 is for 3 phase processing (for casting, filling and modelling).

TECHNOBASE 8000 processes 1-phase and can be modelled directly after mixing.

When it comes to modelling TECHNOBASE does not become sticky and can be modelled to any form desired without the use of aids. The curing time is about 4 – 6 min. according to mixing ratio and product temperature. Once cured TECHNOBASE 6000 has a non-sticky surface. TECHNOBASE 8000 is resistant to acids, alkalis and external impurities and can be dissolved in acetone. In a fluid and doughy state it attacks rubber. Polyethylene and Hostaphan films are not attacked, and can therefore be used to insulate and cover the dough. After completion of polymerisation TECHNOBASE 6000/8000 becomes hard and abrasion-resistant allowing for mechanic machining (grinding, sawing, drilling, milling etc.).

Processing

Preparation

The surface must be roughened (for example, using a hoof rasp, hoof knife or abrasive disc).

The parts getting in touch with TECHNOBASE 6000/8000 (claw, wooden block) must be free of dust and grease and completely dry (we recommend e.g. spirit or alcohol). To achieve reasonable curing times of approx. 4 – 6 min., it may prove practical to cool the mixing vessel and material in summer or to heat it (to about 20 – 23 °C) in winter. At very low outside/product temperatures it is recommended adding TECHNOBASE AKTIVATOR to the mixture.

Processing with ACTIVATOR

TECHNOBASE AKTIVATOR is a co-initiator mixture developed to shorten polymerisation of TECHNOBASE 6000/8000 at low outside/product temperatures.

When outside temperatures are over +10 °C mix TECHNOBASE 6000/8000 liquid and TECHNOBASE AKTIVATOR in a ratio of 14 : 1. At temperatures below +10 °C the mixing ratio is 7 : 1. This corresponds to 500 ml of TECHNOBASE 6000/8000 and 35 ml respectively 70 ml of TECHNOBASE AKTIVATOR.

Mixing TECHNOBASE 6000 – Mixing ratio of powder/liquid about 2 : 1

Shortly before use pour liquid in a mixing beaker (see graduation). Add powder until the mass reaches a still pourable consistency. The mixing ratio can be varied according to individual requirements: Add more powder to obtain a more solid consistency and more fluid for a thinner mixture. Directly after mixing TECHNOBASE 6000 is still fluid enough to be poured. Then the resin turns into a dough-like consistency and becomes ready for filling, and a bit later, for modelling.

Mixing TECHNOBASE 8000 – Mixing ratio of powder/liquid about 3 : 1

Shortly before use pour fluid in a mixing beaker (see graduation). Quickly stir in powder until a non-sticky, moldable paste will be obtained. When dosing individually make sure not to add too much powder. Otherwise the paste will become too dry and cannot be modelled properly.

Application

Spread a thick layer of TECHNOBASE 6000/8000 on the claw sole and wood block. Having pressed the wood block hard onto the claw sole the surplus resin has to be modelled over the sides of the claw and the wood block. The circumferential groove of the wood block has to be filled properly with resin to provide good stability. After about 4 – 6 minutes according to mixing ratio and product temperature TECHNOBASE 6000/8000 will cure completely and once cured will withstand high loads.

Curing

The curing process begins when powder and liquid are mixed. Curing time can be shortened by heating or prolonged by cooling, but not interrupted. Curing brings forth heat. The resulting polymerisation temperature depends on the quantity of TECHNOBASE 6000/8000 added. To avoid possible tissue damages, the resin can be sprinkled with cold water until curing is completed. The synthetic resin can be cut with a hot spatula (caution – flammable!). Openings produced in this way can be easily and permanently closed with freshly-mixed resin. When the healing process has been completed, the plastic will be mechanically removed with the use of appropriate tools (pliers, hammer).

Application

Claw treatment:

Claws are being immobilized during inflammatory diseases by fitting a wood block to the unaffected neighbouring claw using TECHNOBASE 6000/8000 as connector.

Bone Treatment:

To fix the fractured ends in their original position until healing is completed, the bone ends are connected by percutaneous bone screws (Becker) or Kirschner drill wires and individually shaped bridge splints made of TECHNOBASE 6000/8000. Advantages: The low weight is favorable for small domestic animals and furthermore, there is the possibility of repositioning the fragments as long as the callous hasn’t hardened yet (by separating and renewed bonding over the length necessary to ensure bone healing).

Danger/Safety advices

Please follow the instructions on the product packaging and safety data sheets.

Storage advice

Do not store the material above 25 °C. Store in a tightly sealed container. Do not use the material after the expiry date. Keep out of the reach of children.

Field of responsibility

Our usage instructions in spoken and written form, resulting from test carried out by us, are made to the best of our knowledge. However, they remain without obligation, also with reference to any trade mark rights of a third party. This does not free you from checking the products delivered by yourself that they are fit for the intended application and use. We accept no liability with regard to the application, use and processing of the product, since such use is beyond our control remaining exclusively in your field of responsibility. Of course we guarantee for the excellent quality of our products in conformity to our General Terms and Conditions of Delivery.

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