M III®

**Verdünnern für Ebersamen
Culture Medium for Boar Semen
Diluyente de Semen Porcino
Культуральная среда для спермы хряков**

DE

EN

ES

RU

13515/xxxx

31351/5000

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M III® is a culture medium in which boar semen can be preserved at +17°C for a medium storage period. It is delivered in powder form; for the preparation of the final extender sterile pure water is required.

Note: Keep out of the reach of children. Not suitable for the consumption or the treatment of human beings or animals. Do not use for injections. Do not touch the powder without gloves and do not inhale it. Store below +15°C in dark and dry conditions.

1. Composition

M III® consists of the following components:

- Glucose
- Sodium citrate
- EDTA
- Sodium bicarbonate
- Antibiotics

2. Package sizes and alternatives

M III® with Gentamicin sulphate:

- 13515/0001 for 1 litre
- 13515/0005 for 5 litres
- 13515/0010 for 10 litres
- 13515/0083 for 83 litres

100 ml of the prepared extender contain 0.03 g Gentamicin sulphate

3. Application

3.1. Preparation of the final extender

Heat the sterile pure water to a temperature of +32°C to +38°C.

Dissolve 60 g of the M III® powder in 1000 ml of sterile pure water by agitating or stirring.

Allow 20 minutes adaptation time for the pH value to stabilise before starting the dilution of the ejaculate.

3.2. Water quality

The quality of the water used in all the procedures is of a great importance. The water must be purified either by distillation, deionisation, reverse osmosis, or preferably a combination of these methods, and if necessary by additional sterilization. The purified water must have the following characteristics:

Conductivity less than 5 µS/cm

Bacteria level below 1 CFU (colony forming units) per 10 ml

Of utmost importance, all potential sources for bacterial contamination of the prepared water must be identified and eliminated. Unsterile tubing, pipes, or storage containers represent some potential sources for contamination.

If there is no appropriate equipment for water treatment available, bi-distilled and sterile water should be bought from a reliable source.

Never use distilled water for household use!

3.3. Recommended application

The freshly collected boar semen is maintained at a temperature of +30°C to +33°C and analysed. Volume, concentration and motility are determined. Depending on the number of progressive motile sperm cells in the ejaculate, the required amount of extender is calculated (minimum dilution rate 1+4, in order to reach minimal inhibition concentration of antibiotics in the diluted semen). It is recommended to use a diluted sample of the ejaculate for the motility analysis. The extender and the semen must have the same temperature at the moment of the first dilution (+/-1°C).

After the dilution, the semen is no longer maintained in warm conditions. Filling into tubes or bottles of 85 to 100 ml takes place at ambient temperature. After cooling down to ambient temperature, the insemination doses are kept at +17°C.

A daily quality analysis until the application is very much recommended.

In this case semen samples must be pre-warmed for 20 minutes to +37°C prior to the analysis, as the max. motility only is reached after this adaptation phase.

4. Storage, shelf life and important hints

4.1. Storage, shelf life

The airproof sealed package with the powder mixture must be stored in cool (below +15°C), dry and dark conditions. The minimum shelf life is printed on each package.

Ready-made extender can be kept up to 2 days if stored in the refrigerator (+5°C).

4.2. Trouble shooting

A couple of factors are crucial for good semen conservation and insemination success, among others:

- Usage of exclusively sperm friendly materials and surfaces for semen and extender (collection gloves, semen and extender bags, glassware free of mineral residues).
- Good temperature management, particularly avoiding temperatures below +30°C for undiluted semen and below +15°C for diluted semen, as well as temperature fluctuations.
- Avoiding bacterial contamination (appropriate management of boars and preparations of collection; semen collection in clean ambient and with double glove, appropriate hygiene protocol in the lab).
- Avoiding sperm toxic factors (water, UV light, dust, detergent).
- Quality of the purified water (see point 3.2).

4.3. Consulting and training

You can rely on the Minitube scientific experts for advice and technical support concerning all aspects of your semen production and/or artificial insemination process. Workshops and seminars are offered regularly and also on an individual basis. Please refer to our web site for information regarding our current program.

5. Production standard and quality control

5.1. Raw material quality

All raw materials used in Minitube media are produced according to GMP and DIN ISO 9001:2015 norms, and are certified for meeting Ph Eur, BP or USP standards. They are tested and comply with the current international guidelines regulating the quality of pharmaceutical substances.

5.2. Production standard

Production of all Minitube media is performed under GMP conditions in a clean room and documented with weighing records. Under expert veterinary supervision, each batch of the final product is subjected to a series of chemical and physical tests and used under practical conditions to conserve semen. This ensures that Minitube media consistently provide semen conservation with high spermatological efficacy over a wide range of climatic and process conditions.

5.3. Certificates

The general quality certificate and the batch certificate are available upon request.

6. Literature

JOHNSON L.A., WEITZE K.F., FISER P., MAXWELL W.M.C., 2000. Storage of boar semen. Anim. Reprod. Sci. 62, 143 –172.

R. Frangež, T. Gider, M. Kosec, 2004. Frequency of Boar Ejaculate Collection and its Influence on Semen Quality, University of Ljubljana, Veterinary Faculty, Ljubljana, Slovenia.

MINITÜB Protocol, 2007. How to process boar semen on farm. Available upon request.

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All rights reserved, including the translation.
Reprint, even in extracts only with a prior written permission. No part of this manual may be reproduced or transmitted in any form or by any means electronic, mechanical, photocopying, without the prior written permission of Minitüb GmbH.
All brand product names are trademarks or registered trademarks of the respective title holder.
Errors and technical alterations excepted.
Use product only in perfect condition and in compliance with the manual.
Keep safe for future consulting (according to EN ISO 12100).
Translated from German.

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